
Vitapep

Dispositivo per la ginnastica respiratoria
Device for respiratory gymnastic

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Symbols

	Warning: read carefully the warnings contained in the instructions for use
	Manufacturer
	Latex free
	Phthalates free
	Warning: read carefully the instructions for use
	Temperature limits from -25°C to 50°C
	Relative humidity limits from 30% to 95%
	Atmospheric pressure limits from 760 hPa to 1050 hPa
	Production date
	Reference number
	Lot number
	Medical device

Warnings

- Read carefully the instructions for use and keep them in a safe place.
- There are no contraindications to the use of the device
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established
- The use of the device by a child must always take place under the supervision of an adult who is aware of these instructions.
- The device is personal, for one patient only. The device can be used on the same patient for an indefinite number of treatments. Check after every treatment if the device is complete and not damaged, if the device shows signs of damaging or degradation replace it with a new one.
- The manufacturer has defined a useful life of at least 5 years from the first use.

Functioning principle

Among the single alveoli, the bronchi and the distal portions of the lung there are collateral canals. Such canals are open under normal conditions, particularly during deep inhaling phases such as sighing which

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- whilst no resistance must be observed during the inspiratory phase, this must be observed instead during the expiratory phase.

NB: some accessories may not have been included in this manual.

Use of the device

Treatment mode. The device, which must be assembled as described above, must be applied to the patient and the mask must adhere perfectly to the face in order to avoid losses of exhaled air between the mask and the face itself. If the mask dimensions are not suitable for the face of the patient, the device (or at least the mask) must be replaced with a suitable one. The treatment consists of deep inspirations and subsequent expirations, after 10-15 breathing cycles, the device must be disconnected from the patient's face who must alternatively cough and treat out in a forced manner to expel the secretions and thus free the lungs. The gauged hole in the specific patient producing during the expiratory phase the necessary positive pressure for treatment: efficacy.

This is carried out by the doctor using a specific Vitapep device equipped with a pressure measuring instrument. For information purposes, the necessary pressure for treatment efficacy is as follows:

- in common cases: 10-15 hPa
- in emphysema cases: 5-8 hPa
- in atelectasis cases: 12-20 hPa

Note: 1 hPa is equivalent to 1 cm H₂O

Treatment frequency

For treatment frequency, reference to the medical prescription must be made. For information purposes, the recommended frequency is as follows:

- in chronic cases: 10-15 minutes 3/4 times per day;
- in acute illness cases: more often, according to secretion quantities;
- in atelectasis and pharyngitis: every hour

Duration of the treatment

Treatment duration must also be prescribed by the doctor. For chronic patients, the treatment may last for the entire life, but during the periods in which the secretions problems are not observed, the treatment may be stopped. During post-operative treatments, the treatment must be continued through the entire period in which the patient is immobile.

Patient posture

The best position from an efficacy point of view of the treatment is when the patient is seated on a chair with his elbows on the table and keeping his torso erect. For bed-ridden patients, use a slanting back stand with arm rest for the elbows.

Note: the treatment will not produce the

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expected efficacy is:

- the mask does not adhere perfectly to the face;
- expiration is too forced as this would reduce lungs volume;
- frequency and duration do not follow the medical prescription.

Cleaning and maintenance

After each use, cleaning of the components: any non-corrosive detergent product may be used. Components cleaning has to be performed after every use. Vitapep and its component must not have hot or autoclave sterilization.

we recommend dismounting and cleaning periodically (at least once a week) the disk with the gauged holes taking great care in positioning the "O" gasket correctly in their proper housing (see figure 4); follow the instructions described above for subsequent selection of the prescribed hole.

Disposal

All components must be disposed in compliance with current regulations in plastic containers. The table shows the plastic materials used. All the materials used do not contain phthalates and the presence of natural latex is not detected.

COMPONENT	MATERIAL
Body with rotating disk support	ABS
Rotating disk	ABS
Fixing screw	Acciaio
Connector with unidirectional valves	ABS and silicone
Fixing OR gasket	Silicone
Key	ABS and steel

Working and storage conditions

- Temperature from -25°C to +50°C
- Relative humidity from 30% to 95%
- Atmospheric pressure from 760 hPa to 1050 hPa

Air Liquide Medical Systems S.r.l.

All rights reserved. All information contained in this publication may not be used for purposes other than the original ones. Partial or entire reproduction of this publication is forbidden without the written consent of Air Liquide Medical Systems S.r.l. With a view of maintaining and improving the standards of production, functions and increased reliability, the apparatuses of Air Liquide Medical Systems S.r.l. are periodically revised and updated. For this reason, the content of this publication are subject to modification without prior notification. It is possible to obtain the updated version of this document should it become obsolete, by contacting the local distributor or the area agent and specifying the edition in your possession.