EN Light foldable wheelchair
F Fauteuil roulant léger pliant
ES Silla de ruedas ligera plegable
P Cadeira de rodas leve dobrável

Single cross painted aluminium frame

Series **CP770**

Helios Smart



Optional accessories



Cod. **CPA775**
Height adjustable armrest (height from ground 73-81 cm) - PAIR



Cod. **CPA771**
Aluminium height adjustable footrests



Cod. **CPA501-2**
Linear table
55x32 cm (except size 50)
Cod. **CPA504**
Linear table for lifting armrest



Cod. **CPA506-2**
Table for wheelchair
58x38 cm (except size 50)
Cod. **CPA510**
Shaped table for CPA775 lifting armrest
58x38 cm (except 50 cm size)



Cod. **CPA520**
Table for wheelchair
68x57cm (only size 50)



Cod. **CPA777**
Safety wheels



Cod. **CPA772**
Narrow passages/lift wheels

Spare parts

Cod. **CPR774** PAIR
Complete armrest

Cod. **CPR780** PAIR
Detachable footrest

Cod. **CPR241** SINGLE
Plastic footrest only

Cod. **CPR769** PAIR
PU Back wheel Ø 60 cm

Cod. **CPR772** SINGLE
PVC front wheel - ø 20 x 2,5 cm

Cod. **CPR778** PAIR
Brake

Cod. **CPR106** PAIR
Only PU armrest padding

Technical features

SPECIFICATIONS	
Front wheels	PVC, black, steel fork
Rear wheels	Aluminium radius, PU black wheels, aluminium brake
Upholstery	Padded Nylon – Fire reaction tested, according to EN1021-1/2

WEIGHT AND DIMENSION

Measures in cm																	Max load. weight
CP770-40								40		59							
CP770-43								43		62							
CP770-46	106	81	41	90,5	73	51	40÷52	46	44	64	32	15,3	60	20	optional	optional	120
CP770-48								48		67							
CP770-50								50		69							

ARDEA one