EN Cushions
F Coussins
ES Cojines
P Almofadas

opera

With removable cover

Cod. **ST332-44**
Round ø 44 cm

Cod. **ST335-41**
Oval - W: 41 cm

Cod. **ST335-48**
Oval - W: 48 cm









Indicated to give relief in case of post-operative care, decubitus sores, hemorrhoids and prostatitis. Usable on any flat surface and suitable for long car journeys. Hand washable cover at 30°C.

Technical features

Product codes	ST332-44	ST335-41	ST335-48
SPECIFICATIONS			
Inner	moulded polyurethane		
Cover	100% polyurethane		
DIMENSIONS			
	W:	D:	H:	∅
ST332-44	-	-	-	44 cm
ST335-41	41 cm	32 cm	8 cm	-
ST335-48	48 cm	38 cm		-

Inflatable

Cod. **ST304** *ø 40 cm*

Cod. **ST305** *ø 42,5 cm*

Cod. **ST306** *ø 45 cm*



Indicated to give relief in case of post-operative care, decubitus sores, hemorrhoids and prostatitis. Usable on any flat surface and suitable for long car journeys.

Technical features

Product codes	ST304	ST305	ST306
SPECIFICATIONS			
Material	rubber		
DIMENSIONS			
∅	40 cm	42,5 cm	45 cm